



## The DISCO bike keeps you independent and mobile, giving you the freedom to exercise and a means of transport for doing daily errands when it suits you

DISCO is recommended for people with special needs or who are uncomfortable mounting/dismounting a normal bike, as well as anyone who just wants a different kind of tricycle. The comfortable sofa chair at low height guarantees a comfortable and safe bike trip, every time!

The model is available as a disability tricycle for both adults and children. The DISCO tricycle is renowned for comfort, durability and excellent performance:

- easy to mount and safe to ride
- comfortable chair seat which can be easily moved backwards or forwards while riding
- extremely rugged with a long lifetime (10+ years)
- 7 gears with or without reverse function
- adapted to your special needs and delivered within 2-3 weeks

DISCO models are also available with an electric motor, and all bikes can be adapted to your special needs.

### Standard equipment:

- 7 speed Shimano Nexus
- Differential
- Foot brake
- Hand disc brakes front and back
- Parking brake
- Adjustable seat and handlebar
- Front basket
- Lock w/ 2 keys
- Rear view mirror
- Seat cushions
- Tyres:

Disco: Schwalbe Pick-Up

Disco Medi: Schwalbe

Marathon Plus

# DISCO P25:

(front wheel motor)

The Disco P25 has a speed sensor, which ensures that the motor only operates when you are pedaling. You can choose between 5 levels with respect to how much assistance you require from the motor. For help when starting to pedal, you might want to make use of the built-in walk assist (start help) by turning the throttle (up to max 6 km/h). You can drive this bike alone by means of the assist levels, but you can also use the throttle to “override” the current assist level and give full power e.g. when you climb hills or need to pass someone or something quickly. When you let go of the throttle, the bike falls back into the chosen assist level. The motor automatically cuts out at max 25 km/h. With a fully charged standard battery, the bike has a range of up to 35 km. Depending on how much pedaling you do, chosen assist level as well as wind and road conditions, the range can reach distances well above that level.

The display on the handlebar shows you information such as speed, distance, battery level and assist level. The battery is removable, so charging can be done where you prefer.

### DISCO P25 gear options:

- 7 speed free wheel
- 7 speed with back pedal and reverse gear



DISCO P25



DISCO Medi P25 (16" wheels)



DISCO P25

DATA:	Disco Medi	Disco	Disco colour options:	
Length:	171 cm	185 cm		Matte red
Width:	68 cm	75 cm		Matte blue
User measure inside leg	65 cm	72 cm		Matte green
User weight, max:	125 kg	125 kg		Matte grey
Mounting height:	19.5 cm	24.5 cm	<i>PF mobility bikes are painted using an environment-friendly powder lacquer</i>	
Wheel size:	16"	20"		
Motor:	250W	250W		
Battery:	396 Wh or	692 Wh		